

Vertical Pressure Steam Sterilizer (Digital Display Automation)



Global B2B Marketplace Focused on the Supply and Procurement of Medical Equipment and Devices

w w w . m e d w i s h . c o n

Technical Parameter

TECHNICAL DATA

Fully stainless high-alloy steel structure

Digital display of working status, touch of key

Auto discharge the cool air, and steam discharging automatically after sterilization

Automatically shut off with beep reminding after sterilization

With two stainless steel sterilizing baskets

Easy to operate, safe and reliable

Model:	LS-75LD
Chamber volume	75L(φ400×600)mm
Working pressure	0.22MPa
Working temperature	134°C
Max. working pressure	0.23 Mpa
Heat average	≤±1°C
Timer scope	$0\sim99$ min $/0\sim99$ hour 59 min
Adjustment of temperature	105~134°C
Power	3*1.5KW /AC220V.50Hz
Overall dimension	550×530×990(mm)
Transport dimension	650×650×1150(mm)
G.W/N.W	90Kg/67Kg











Follow Us









Jiangsu Medwish Network Technology Co., Ltd